### INSTRUCTIONS FOR USE

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## EN

### INSTRUCTIONS FOR USE

HandX Needle Holder instrument IFU  
HandX Fenestrated Grasper instrument IFU

4 Meir Ariel St., Grand Netter Building Netanya 4250574, Israel Telephone +(972)-77-3630300

Obelis S.A Registered Address: Bd. Général Wahis, 53 1030 Brussels, Belgium

0123

### HandX DEVICE DESCRIPTION

The Human Xtensions HandX is a handheld powered laparoscopic device.

The HandX is electromechanically controlled, includes hardware and software and comprises two major parts, the Handpiece and the instrument (Figure 1).

### Main HandX Components as Illustrated in Figure 1

- 1. Arc & Spacer
- 2. CI (Control Interface)
- 3. Finger Unit
- 4. Handpiece Body
- 5. Gearbox
- 6. Gearbox Cap
- 7. Shaft
- 8. Articulating Tip
- 9. 5mm Needle Holder
- 10. 5mm Fenestrated Grasper
- 11. Finger Pads
- 12. Power Adapter
- 13. Sterile Power Cable
- 14. Protective Plug
- 15. Release Key

### INSTRUMENT

- Articulating Tip
- Shaft
- Gearbox Cap
- Protective Plug
- Gearbox
- Release Key

### INDICATION

The Human Xtensions HandX™ device is indicated for surgical treatment of tissue and organ inflammation, benign and/or malignant space occupying lesions, and anatomic deformations in the abdominal, thoracic, and pelvic cavities.

### INTENDED USE

The Human Xtensions HandX™ device is intended to facilitate repair and/or removal of organs and soft tissues during laparoscopic procedures through tissue manipulation and mobilization, ie, stitching and grasping using the Human Xtensions HandX™ device.

### PATIENT POPULATION

The target patient population is restricted to paediatrics and adults suitable for laparoscopic surgery.

### CONTRAINDICATIONS

Where laparoscopic surgery is contraindicated.

### WARNINGS

1. The disposable components of the device are provided STERILE and are intended for use in a single procedure only. DISCARD AFTER USE DO NOT RESTERILIZE. Reuse reprocessing, or re-sterilization may compromise the structural integrity of the device and/or lead to device failure that in return may result in patient injury, illness or death.
2. Endoscopic procedures should be performed only by physicians having adequate training and familiarity with endoscopic techniques. A thorough understanding of the operating principles, risks versus benefits, and the hazards involved in an endoscopic approach is necessary to avoid possible injury to the user and/or patient.
3. HandX™ instrument cannot be used by itself and are intended to be used only with a HandX™ Handpiece.
4. Prior to use, visually inspect the packaging/pouch/blister for breaches of packaging integrity. Do not use the HandX Handpiece, instruments, or HandX Disposables if the device or package is damaged/opened/dropped, as it could compromise the safety, functionality, and sterile barrier of the HandX system.
5. Do not use past the expiry date.
6. Used HandX instruments should be disposed according to bio-hazard disposal practices in accordance with the facility disposal procedure.
7. Do not connect wet instruments to the HandX Handpiece.
8. No modification of this device is allowed without authorization of the manufacturer.

NOTE: Defibrillation was tested and in compliance with IEC 60601-1.

### USING THE DEVICE

Follow Pre-operational assembly per HandX™ Instructions for Use.

### TRANSPORT CONDITIONS

- HandX instruments:
- Temperature: -10°C to 50°C
  - Humidity: 10% - 85%

### STORAGE CONDITIONS

- HandX instruments:
- Temperature: 15°C to 30°C
  - Humidity: 10% - 85%

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