EN

FR

IT

DE

ES

GR

INSTRUCTIONS FOR USE

HandX Needle Holder instrument IFU

HandX Fenestrated Grasper instrument IFU

HandX Self Righting Needle Holder instrument IFU

This page is intentionally left blank.

Figure 1: Main HandX Components

Do not use if package is damaged or open | Ne pas utiliser si l'emballage est endommagé ou ouvert | Non utilizzare se la confezione è danneggiata o aperta | Inhalt bei beschädigter oder offener Verpackung nicht verwenden | No usar si el envase está dañado o abierto | Να μη χρησιμοποιείται αν η συσκευασία έχει υποστεί φθορά ή έχει ανοιχτεί

Follow instructions for use | Consulter le mode d'emploi | Seguire le Istruzioni per l'uso | Befolgen Sie die Gebrauchsanweisung | Seguir las instrucciones de uso | Ακολουθήστε τις οδηγίες χρήσης

Defibrillation-Proof Type BF Applied Part | Pièce en contact avec le patient de type BF protégée contre la défibrillation | Parte applicata di tipo BF a prova di defibrillazione | Anwendungsteil vom Typ BF, defibrillations-geschützt | Componentes aplicados a prueba de desfibrilación- de tipo BF | Εφαρμοζόμενο μέρος ανθεκτικό στην απινίδωση τύπου BF

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HandX Needle Holder instrument IFU

HandX Fenestrated Grasper instrument IFU

HandX Self Righting Needle Holder instrument IFU

HumanTouch Surgical Ltd.
Grand Netter Building,
4 Meir Ariel st, Netanya, 4250574. Israel
Telephone: +(972) 77 3630300

MedNet EC-REP IIb GmbH, Borkstrasse 10,
48163 Münster, Germany

HandX DEVICE DESCRIPTION

The Human Xtensions HandX is a handheld powered laparoscopic device.

The HandX is electromechanically controlled, includes hardware and software and comprises two major parts, the Handpiece and the instrument (Figure 1).

Main HandX Components as Illustrated in Figure 1

- | | | |
|---------------------------|-------------------------------------|-------------------------|
| 1. Arc & Spacer | 7. Shaft | 12. Finger Pads |
| 2. CI (Control Interface) | 8. Articulating Tip | 13. Power Adapter |
| 3. Finger Unit | 9. 5mm Needle Holder | 14. Sterile Power Cable |
| 4. Handpiece Body | 10. 5mm Fenestrated Grasper | 15. Protective Plug |
| 5. Gearbox | 11. 5mm Self-Righting Needle Holder | 16. Release Key |
| 6. Gearbox Cap | | |

INSTRUMENT

- | | | |
|--------------------|-------------------|---------------|
| • Articulating Tip | • Gearbox Cap | • Gearbox |
| • Shaft | • Protective Plug | • Release Key |

INDICATION

The Human Xtensions HandX™ device is indicated for surgical treatment of tissue and organ inflammation, benign and/or malignant space occupying lesions, and anatomic deformations in the abdominal, thoracic, and pelvic cavities.

INTENDED USE

The Human Xtensions HandX™ device is intended to facilitate repair and/or removal of organs and soft tissues during laparoscopic procedures through tissue manipulation and mobilization, ie, stitching and grasping using the Human Xtensions HandX™ device.

PATIENT POPULATION

The target patient population is restricted to paediatrics and adults suitable for laparoscopic surgery.

CONTRAINDICATIONS

Where laparoscopic surgery is contraindicated.

WARNINGS

1. The disposable components of the device are provided STERILE and are intended for use in a single procedure only. DISCARD AFTER USE DO NOT RESTERILIZE. Reuse reprocessing, or re-sterilization may compromise the structural integrity of the device and/or lead to device failure that in return may result in patient injury, illness or death.
 2. Endoscopic procedures should be performed only by physicians having adequate training and familiarity with endoscopic techniques. A thorough understanding of the operating principles, risks versus benefits, and the hazards involved in an endoscopic approach is necessary to avoid possible injury to the user and/or patient.
 3. HandX™ instrument cannot be used by itself and are intended to be used only with a HandX™ Handpiece.
 4. Prior to use, visually inspect the packaging/pouch/blister for breaches of packaging integrity. Do not use the HandX Handpiece, instruments, or HandX Disposables if the device or package is damaged/opened/dropped, as it could compromise the safety, functionality, and sterile barrier of the HandX system.
 5. Do not use past the expiry date.
 6. Used HandX instruments should be disposed according to bio-hazard disposal practices in accordance with the facility disposal procedure.
 7. Do not connect wet instruments to the HandX Handpiece.
 8. No modification of this device is allowed without authorization of the manufacturer.
 9. Do not connect wet instruments or HandX Disposables to the HandX Handpiece.
 10. The HandX Handpiece shall maintain control of the articulation tip movement at all times. In case of an irregular articulation tip movement, stop using the device.
 11. Do not disconnect the instrument during operation. In case of inadvertent instrument disconnection, replace the instrument.
 12. In addition to the listed warnings, the use of the HandX system may cause/lead to the following side effects:
 - a Uncommon side-effects (1-10 out of 1,000):
 - i Inability to perform the intended surgical task due to a damaged or non-functional device.
 - ii Bleeding or tissue damage due to un-intentional HF energy burst.
 - iii Damage to the trocar seal, loss of intra-abdomen pressure during surgery.
 - b Rare side-effects (1-10 out of 10,000):
 - i Patient inflammation or infection caused by a reaction from foreign materials.
 - ii Bleeding and tissue damage due to loss of performance or tip control.
-  Potential adverse events- Any serious incident must be reported to the manufacturer and the competent authority of the member state.

NOTE: Defibrillation was tested and in compliance with IEC 60601-1.

USING THE DEVICE

Follow Pre-operational assembly per HandX™ Instructions for Use.

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